


LOHMANN

X-RAY

80/6 EP



80 KV X-RAY INSERT
OPEN ANODE
DIRECTIONAL BEAM

07094

tube voltage max.:	80 kV	nom. focal spot:	1,8 x 1,8 mm
tube current max.:	6 mA	inherent filtration:	1,5 mm glass
(depends on cooling)		anode angle:	25°
tube power max.:	480 Watt	beam angle:	80° x 60°
filament voltage ca.:	4,0 V	cooling:	oil, gas
filament current:	2,8 ± 0,2 A	frequency:	50 - 60 Hz

Version
84008

11

German Craft and Technology since 1946